

Western 转移缓冲液(NC,PVDF)

产品简介:

Western 转移缓冲液又称蛋白印迹转移缓冲液或转膜缓冲液,适用于 Western Blot 实验中把蛋白质转移到 PVDF 膜或硝酸纤维素(NC)膜上,以便进行抗体的孵育和下游实验。

Leagene Western 转移缓冲液(NC,PVDF)主要成分为 Tris、甘氨酸、甲醇等,不适用于把蛋白质转移到尼龙膜上,该试剂为 1×溶液,直接使用。该试剂仅用于科研领域,不适用于临床诊断或其他用途。

产品组成:

名称	编号	PW0034	Storage
Western 转移缓冲液(NC,PVDF)		500ml	RT
使用说明书		1 份	

操作步骤(仅供参考):

- 1、根据实验具体要求操作。
- 2、大多数蛋白质转移,应置于 4℃进行,以防止蛋白降解。

注意事项:

- 1、密闭保存,开启后尽快用完。
- 2、为了您的安全和健康,请穿实验服并戴一次性手套操作。
- 3、试剂开封后请尽快使用,以防影响后续实验效果。

有效期: 12 个月有效。

相关产品:

产品编号	产品名称
DC0032	Masson 三色染色液
PS0013	RIPA 裂解液(强)
PE0103	Acr-Bis(30%,29:1)
PW0053	Western 抗体洗脱液(碱性)
PW0060	通用显影液
PW0084	丽春红 S 染色储存液(10×Ponceau S)
TO1013	丙二醛(MDA)检测试剂盒(TBA 比色法)

文献引用:

- 1、Wang Yanbo,Ye Xiaojun Su Wenwen,et al.Diosmin ameliorates inflammation, apoptosis and activates PI3 K/AKT pathway in Alzheimer' s disease rats.METABOLIC BRAIN DISEASE.August 2024.10.1007/s11011-024-01388-7.(IF 3.2)

- 2、 Hua Zhen,Dai Shijie,Li Shaoshuo,et al.Deciphering the protective effect of Buzhong Yiqi Decoction on osteoporotic fracture through network pharmacology and experimental validation.Journal of Orthopaedic Surgery and Research.February 2023.10.1186/s13018-023-03545-7.(IF 2.6)

注：更多使用本产品的文献请参考产品网页