

碱性琼脂糖凝胶电泳缓冲液(10×)

产品简介:

碱性琼脂糖凝胶电泳缓冲液(10×)是一种 10 倍浓缩的碱性 DNA 电泳缓冲液, 由 NaOH、EDTA 等组成, 呈强碱性, 使用时稀释至 1×使用, 主要用于碱性琼脂糖凝胶电泳。该试剂仅用于科研领域, 不适用于临床诊断或其他用途。

产品组成:

名称 \ 编号	NA0070	Storage
碱性琼脂糖凝胶电泳缓冲液(10×)	500ml	RT
使用说明书	1 份	

操作步骤(仅供参考):

- 1、用去离子水稀释碱性琼脂糖凝胶电泳缓冲液(10×)至 1×。
- 2、用 1×碱性琼脂糖凝胶电泳缓冲液制备用于碱性琼脂糖凝胶。
- 3、加样品于琼脂糖加样孔, 电泳。

注意事项:

- 1、如果每次的使用量很小, 可以适当分装后再使用。
- 2、为了您的安全和健康, 请穿实验服并戴一次性手套操作。
- 3、试剂开封后请尽快使用, 以防影响后续实验效果。

有效期: 12 个月有效。

相关产品:

产品编号	产品名称
DH0006	苏木素伊红(HE)染色液(醇溶)
NA0030	Tris-乙酸电泳缓冲液(50×TAE)
NH0043	SSC 缓冲液(20×, pH7.0)
NR0040	RNase A(10mg/ml)
OR0001	pH 标准缓冲溶液(pH=4.00)
PT0001	BCA 蛋白定量试剂盒
PW0040	Western blot 一抗稀释液
TC1243	甘油三脂(TG)检测试剂盒(GPO-PAP 单试剂比色法)

文献引用:

- 1、 Yaping Gao,Yun He,Liyi Chen,et al.Chimeric Phi29 DNA polymerase with helix–hairpin–helix motifs shows enhanced salt tolerance and replication performance.Microbial Biotechnology.May 2021.10.1111/1751-7915.13830.(IF 5.813)
- 2、 Yaping Gao,Yun He,Igor Ivanov,et al.Expression and functional study of VpV262 Pol, a moderately halophilic DNA polymerase from the Vibrio parahaemolyticus phage VpV262.ENZYME AND MICROBIAL TECHNOLOGY.June 2020.10.1016/j.enzymictec.2020.109588. (IF 3.448)

注：更多使用本产品的文献请参考产品网页