

RNase A(10mg/ml)

产品简介:

Ribonuclease A 简称 RNase A, 中文名为核糖核酸酶 A, Leagene RNase A(10mg/ml) 主要由 RNase A、Tris-HCl 等组成, 含有极少量 DNase, 用于消化 RNA, 可用于质粒和基因组抽提等相关实验中去除 RNA, 常用于各种分子生物学、细胞生物学等相关实验, 可直接使用。该试剂仅用于科研领域, 不适用于临床诊断或其他用途。

产品组成:

名称	编号			Storage
	NR0040	NR0040		
RNase A(10mg/ml)	1ml	5ml		-20℃ 避光
使用说明书			1 份	

操作步骤(仅供参考):

1、根据具体实验要求操作, 一般工作浓度为 $1 \sim 10 \mu\text{g/ml}$ 。

注意事项:

- 1、应注意避免污染和反复冻融。
- 2、为了您的安全和健康, 请穿实验服并戴一次性手套操作。
- 3、试剂开封后请尽快使用, 以防影响后续实验效果。

有效期: 12 个月有效。低温运输, -20℃保存。

相关产品:

产品编号	产品名称
DP0013	GUS 染色液(即用型)
NR0042	RNase A(10mg/ml,DNase free)
PE0025	SDS-PAGE 蛋白加样缓冲液(5×)
R00738	乙酸钠溶液(3mol/L,pH5.2,无菌)
DH0005	Mayer 苏木素染色液
TC1167	尿素(Urea)检测试剂盒(脲酶波氏比色法)

文献引用:

- 1、Chao Tian, Meng Wang, Xueqi Shi,et al.Discovery of (2-(pyrrolidin-1-yl)thieno[3,2-d]pyrimidin-4-yl)(3,4,5-trimethoxyphenyl)methanone as a novel potent tubulin depolymerizing and vascular disrupting agent.EUROPEAN JOURNAL OF MEDICINAL CHEMISTRY. May 2022. 10.1016/j.ejmech.2022.114466.(IF 7.088)
- 2、Chao Tian, Xuanchen Chen, Zhili Zhang,et al.Design and synthesis of (2-(phenylamino)thieno[3,2-d]pyrimidin-4-yl)(3,4,5-trimethoxyphenyl)methanone analogues as potent anti-tubulin polymerization agents.EUROPEAN JOURNAL OF MEDICINAL CHEMISTRY.September 2019.10.1016/j.ejmech.2019.111679.(IF 4.833)

注: 更多使用本产品的文献请参考产品网页