

牛血清白蛋白溶液(10% BSA)

产品简介:

牛血清白蛋白溶液(10% BSA)是用于包被、蛋白添加剂等多种用途的试剂, 主要由 BSA(组分 V)、去离子水组成。该试剂仅用于科研领域, 不适用于临床诊断或其他用途。

产品组成:

名称	编号	R00912	Storage
牛血清白蛋白溶液(10% BSA)		100ml	4°C
使用说明书		1 份	

操作步骤(仅供参考):

- 1、根据实验具体要求操作。

注意事项:

- 1、小心操作, 尽量避免细菌污染。
- 2、含有 BSA 的溶液颜色容易发黄, 有泡沫, 一般不影响使用。
- 3、为了您的安全和健康, 请穿实验服并戴一次性手套操作。
- 4、试剂开封后请尽快使用, 以防影响后续实验效果。

有效期: 6 个月有效。低温运输, 4°C 保存。

相关产品:

产品编号	产品名称
CC0005	磷酸缓冲盐溶液(1×PBS,无钙镁)
CC0128	胰蛋白酶-EDTA 溶液(0.25%:0.02%,含酚红)
CS0001	ACK 红细胞裂解液(ACK Lysis Buffer)
DC0032	Masson 三色染色液
TC0699	植物总糖和还原糖检测试剂盒(DNS 比色法)
TC1167	尿素(Urea)检测试剂盒(脲酶波氏比色法)

文献引用:

- 1、Di Wu,Tiancheng Gong,Zhongxiang Sun,et al.Dual-crosslinking gelatin-hyaluronic acid methacrylate based biomimetic PDAC desmoplastic niche enhances tumor-associated macrophages recruitment and M2-like p

olarization.INTERNATIONAL JOURNAL OF BIOLOGICAL MACROMOLECULES.April 2024.10.1016/j.ijbiom
ac.2024.131826.(IF 7.7)

- 2、 Haopeng Pan,Shajun Zhu,Tiancheng Gong,et al.Matrix stiffness triggers chemoresistance through elevated autophagy in pancreatic ductal adenocarcinoma.Biomaterials Science.October 2023.10.1039/D3BM00598D.
(IF 6.6)
- 3、 Yue Zhang,Ruijiao Chang,Muqiong Li,et al.Docetaxel-loaded lipid microbubbles combined with ultrasound-t
riggered microbubble destruction for targeted tumor therapy in MHCC-H cells.OncoTargets and Therapy.N
ovember 2022.10.2147/OTT.S102101.(IF 4.345)
- 4、 Feng Ouyang,Li Zhao,Qi Shuai.Photothermally driven degradable nanomissile for single-wavelength photot
herapy and hypoxia-activated chemotherapy.INORGANIC CHEMISTRY COMMUNICATIONS.August 2023.
(IF 3.8)

注：更多使用本产品的文献请参考产品网页